



Canadian Hard of Hearing Association
Association des malentendants canadiens



Get in the **Hearing Loop**
Bienvenue à la **boucle auditive**

For Audiologists and Hearing Aid Providers

Please advise me if my hearing device(s) has a telecoil(s) / T-coil. If they do, please activate and program them for use in a hearing loop. Also, please write instructions for their use on the back of this card. If I don't have a T-coil(s), can my device(s) be retrofitted or is there an accessory remote control / streamer telecoil option?

Program Recommendations

For open RIC fittings: manual T only program

For closed / occluded fittings: manual M + T program in social settings

For closed / occluded fittings: manual T only program for high noise environments like airports or train stations.

Verify mic and telecoil responses **match** in gain, output, and frequency response in user programmed mode (use ANSI-SPLIV test procedure).

I would like to use my T-coil(s) in a hearing loop or with a neckloop and an assistive listening device or system.

Please write clear instructions on how to use my manually accessible telecoil program.

For example: the numbered position of the program, number of beeps or voice prompt, and location of the switch to activate the T-coil.

Please write instructions here:



Audiologist or hearing aid provider name and contact information