



Canadian Hard of Hearing Association
Association des malentendants canadiens



Get in the **Hearing Loop**
Bienvenue à la **boucle auditive**

For Audiologists and Hearing Aid Providers

Please advise me if my hearing device(s) has a telecoil(s) / T-coil. If they do, please activate and program them for use in a hearing loop. Also, please write instructions for their use on the back of this card. If I don't have a T-coil(s), can my device(s) be retrofitted or is there an accessory remote control / streamer telecoil option?

Program Recommendations

For open RIC fittings: manual T only program

For closed / occluded fittings: manual M + T program in social settings

For closed / occluded fittings: manual T only program for high noise environments like airports or train stations.

Verify mic and telecoil responses **match** in gain, output, and frequency response in user programmed mode (use ANSI-SPLIV test procedure).

I would like to use my T-coil(s) in a hearing loop or with a neckloop and an assistive listening device or system.

Please write clear instructions on how to use my manually accessible telecoil program.

For example: the numbered position of the program, number of beeps or voice prompt, and location of the switch to activate the T-coil.

Please write instructions here:



Audiologist or hearing aid provider name and contact information

This resource, its contents and design are used with permission from the Hearing Loss Association of America (HLAA). Copyright 2022. This material is a product of the HLAA Get in the Hearing Loop (GITHL) Program. For more information go to www.hearingloss.org/programs-events/get-hearing-loop/.