

HEAR HERE

This symbol indicates that a **HEARING LOOP** is available.

Activate the telecoil (t-coil) in your hearing device to hear sound with less background noise – directly into your device.

If you do not have a t-coil enabled device, request a headset + receiver.

For questions about hearing loops and hearing loss contact or visit the Canadian Hard of Hearing Association:

loops@chha.ca
www.getinthehearingloop.ca



Hearing Loop
Switch hearing aid to T-coil

A LIFE-CHANGING SOLUTION FOR PEOPLE WITH HEARING LOSS

HEARING LOOPS

- **Simple** - hearing loops reduce background noise and send sound from a microphone directly to t-coil enabled devices
- **Discreet** - there is no need to ask for additional equipment
- **Effortless** - most hearing devices have an external toggle switch to enable the t-coil

Talk to your audiologist about t-coil equipped hearing devices. If you already have one of these devices, ask if it is t-coil enabled and learn how to activate the t-coil program.

This resource, its contents and design are used with permission from the Hearing Loss Association of America (HLAA). Copyright 2022. This material is a product of the HLAA Get in the Hearing Loop (GITHL) Program. For more information go to www.hearingloss.org/programs-events/get-hearing-loop/.



Microphone



Hearing Aids
Cochlear Implants
Bone Conduction

Toggle to access
the t-coil